

# The Gilson Rice Shaker: A Versatile Tool for Specific Gravity Testing



If you perform asphalt theoretical maximum specific gravity (Rice) tests on a routine basis, it's likely you're already familiar with Gilson's Rice Shaker. However, even if you test materials other than asphalt, you should consider if this versatile time saver belongs in your laboratory. The Rice Shaker is the key to accurate, repeatable, and efficient specific gravity determinations for asphalt, aggregates, and soils.

## Specific Gravity/Relative Density: What You Should Know

Specific gravity is a fundamental characteristic of all materials widely used in construction materials testing. The property is applied when proportioning ingredients to mix composite materials like concrete or asphalt for density determinations, bearing capacity and strength of soils in engineered fills, or calculations for aggregates' bulk density, absorption, and void content.

Specific gravity (SG) or the more technically accurate term, [relative density \(RD\)](#), are used interchangeably to describe the ratio of the mass of a material compared to the mass of an equal volume of a reference material, usually water. Specific gravity or relative density directly compares a material's density to water's density, and the resulting ratio is a unitless number.

Our reference substance of water has a specific gravity of 1.0, so materials with greater specific gravities sink, and those with a lower value float. We can bob for apples in water because their relative density is around 0.40, but limestone with an SG of 1.55 sinks like, well, a stone.

## How is Specific Gravity/Relative Density Measured?

The density of water is defined as  $1.0\text{g/cc}^3$  at  $4.0^\circ\text{C}$  ( $39.2^\circ\text{F}$ ). Primarily for convenience, specified water temperatures for SG tests of most construction materials range from  $20^\circ\text{C}$  to  $25.0^\circ\text{C}$  ( $68.0^\circ\text{F}$  to  $77.0^\circ\text{F}$ ), essentially near room temperature. A correction factor for the lower water density is applied when calculating test results. In the simplest test methods, the mass of sample material is determined in air and again when immersed in water:

$$\text{Relative Density} = B/(B-C)$$

B = Sample Mass in Air

C = Sample Mass in Water

Individual specific gravity test methods may include variations in sample preparation, procedure, and calculations to allow further determinations of apparent relative density, average relative density, absorption, and void content.

Large-volume coarse aggregate samples or intact asphalt or concrete specimens are contained in a wire mesh basket or held in a weighing cradle for weighing in water. Baskets or cradles can be suspended from a hook built into the bottom of many scales and balances. The suspended samples are immersed in water for mass and volume displacement determinations. Smaller samples of loose, unbound materials must be contained in flasks or pycnometers to determine the volume and immersed mass.

## Removing Air from Specific Gravity Samples

Air entrapped in the sample can introduce significant errors for any density determination method. For large samples weighed by suspension in water, entrapped air can be avoided by preparing porous materials to a saturated, surface-dry (SSD) condition. Immersed samples contained in pycnometers require greater attention to remove all entrapped air, and some standard test methods require specific practices for deairing by agitation or vacuum application.

The manual deairing process requires gentle, hands-on agitation while the sample is under vacuum to release air effectively. Heating the pycnometer and its contents is sometimes used to promote a boiling action under a partial vacuum. Deairing can be tedious, but neglecting these measures directly impacts the accuracy and repeatability of specific gravity determinations.

## Specific Gravity of Asphalt Mixtures


The Rice test, named after its creator, James Rice, determines the theoretical specific gravity of an asphalt mix without air voids (Gmm). In test methods ASTM D2041 and AASHTO T 209, an asphalt sample is submerged in water. A vacuum is applied to remove free and entrapped air from the container to prepare it for volumetric determination and weighing. Options for the container used in the methods include a metal or plastic pycnometer or a heavy-wall glass vacuum flask.

Many users of the early versions of the Rice test found the reproducibility of the results to be poor. Over time, researchers improved many aspects of the test, including the deairing process, to produce more reliable test results. A device that mechanically agitated the sample for a specified time under a designated vacuum was found to improve test repeatability and, as a bonus, freed technicians from the task of manual agitation.

## Gilson's Innovative Rice Shaker

Early versions of Rice test agitators were crudely adapted abrasive parts cleaners with little or no control over vibration frequency or amplitude and rudimentary on/off switches. Test repeatability was still erratic, and excessive vibration levels easily stripped asphalt from the aggregate and could degrade the aggregate particles.

Gilson's SGA-5R Rice Shaker was conceived as a quality instrument, purpose-built for laboratory specific gravity testing. The mechanical agitator applies and maintains a controlled level of agitation during the vacuum deairing of pycnometers and flasks for a variety of SG tests. The Rice Shaker allows the operator to tend to other tasks during the test cycle.



Gilson builds each Rice Shaker in-house from start to finish. The SGA-5R is powered by a variable-speed electric motor regulated with a rotating variable frequency controller. An eccentric weight attached to the motor causes the base plate to vibrate during operation. Total control of vibration frequency and intensity ensures reliable and repeatable test results, prevents sample damage from abrasion or fractures, and optimizes deairing operations.

A three-position power switch selects timed or manual operating modes. Gilson's exclusive digital countdown timer is built into the front control panel and controls the operation time of the shaker with

±0.25 second accuracy. The four-mode programmable controller features start/stop and pause/resume functions and displays up to 99:59hr: min in hours:minutes, minutes, minutes:seconds, or seconds. A special Texas model, the SGA-5RT, has modified vibration characteristics that meet Texas DOT requirements for test method TEX-227-F.

Gilson's exclusive [EZ-Clamp system](#) secures all of the specific gravity vessels with the push of a button and a quick twist of the clamping knobs. Specific gravity vessels are easily mounted to or removed from the Rice Shaker with lids, fixtures, and vacuum lines in place. An EZ Clamp Upgrade Kit is available to convert older Gilson Rice Shakers to the EZ Clamp system.

## A Device for All Specific Gravity Applications

The asphalt Rice test is not the only procedure in construction materials laboratories that combines the application of a vacuum with agitation for deairing. The same sample immersion and deairing process for weighing and volume determinations is also used in specific gravity tests for aggregates and soil.

[Pycnometers](#) or [filter flasks](#) are positioned on the base plate and secured with the EZ Clamp system on a top plate between two upright rods. The base and top plates accept SG-16A and SG-18A Aluminum Pycnometers, and available adapter sets modify the shaker to fit additional pycnometers and various glass vacuum or volumetric flasks.

NOTE: Review individual test methods to confirm the allowable use of mechanical agitators.

### Pycnometers and Tests Compatible with Gilson Rice Shaker

Pycnometer Type	Gilson Model	Required Adapter	Sample Material	ASTM	AASHTO			
2,000g Aluminum	<a href="#">SG-16A</a>	Included	Asphalt	D2041	T 209			
4,000g Aluminum	<a href="#">SG-18A</a>							
6,000g Plastic	<a href="#">SG-15</a>	<a href="#">SGA-7</a>						
1,000mL, 2,000mL, or 4,000mL Thick-Wall Glass Flasks	<a href="#">GW-74, GW-75, GW-76</a>	<a href="#">SGA-8</a>						
1qt Glass Jar & Top	<a href="#">SG-2</a>					Fine Aggregate	C128	T 84
LeChatelier 250mL Glass Flask	<a href="#">SG-24</a>					Hydraulic Cement	C188	T 133
100mL, 250mL, or 500mL	<a href="#">SG-100, SG-250, SG-500</a>		Soil	D854	T 100			

## Gilson Accessories Improve Precision

Optional accessories monitor vacuum levels in the vessels, and models that regulate vacuum and measure shaker energy are available.

- Gilson's unique [Digital Residual Pressure Manometer](#) adds another layer of accuracy and repeatability. This mercury-free digital manometer displays instant absolute vacuum pressures from 0 to 1,000 mm/Hg with 0.1 mm/Hg readability. The Digital Residual Pressure Manometer is also available with NIST Traceable 3-point Certification.
- [AutoRice™](#) teams up with the Rice Shaker to create the ideal asphalt specific gravity testing system. The digital controller automatically measures and regulates the applied vacuum pressure and time to ensure consistent deairing with minimal attention. The AutoRice™ can also be outfitted with a [vibration sensor](#) to measure the amplitude and frequency of the shaker to establish optimum settings for consistent results.
- The [PumpSaver Air Dryer](#) works with the AutoRice or any vacuum system to eliminate moisture and prevent damage to expensive vacuum pumps.
- [Vacuum pumps](#) or [aspirators](#) aid in deairing samples during agitation and have many other applications around the lab.
- An [Air/Gas Dryer](#) removes moisture from evacuated air. Installing one in the line between the vacuum pump and the pycnometer protects the pump against moisture damage.

We hope this blog post has helped you discover the versatility of our SGA-5R Rice Test Shaker. To discuss your specific application, please contact the [Testing Experts](#) at Gilson.

### Meet the Author



#### *Ben Backus*

*Content Contributor and (retired) Vice President of  
Technical Operations for Gilson Company, Inc.*

Ben is a member of the American Concrete Institute and active on a number of ASTM C09 and ASTM D18 subcommittees. He has over 40 years of experience in the construction materials testing industry and is a NICET senior engineering technician.



### Learn More:

#### Test Standards:

[AASHTO T 84](#), Standard Test Method for Relative Density (Specific Gravity) and Absorption of Fine Aggregate

[AASHTO T 100](#), Standard Method of Test for Specific Gravity of Soils

[AASHTO T 209](#), Standard Method of Test for Theoretical Maximum Specific Gravity and Density of Asphalt Mixtures

[ASTM C128](#), Standard Test Method for Relative Density (Specific Gravity) and Absorption of Fine Aggregate

[ASTM D854](#), Standard Test Methods for Specific Gravity of Soil Solids by the Water Displacement Method

[ASTM D2041](#), Standard Test Method for Theoretical Maximum Specific Gravity and Density of Asphalt Mixtures

[ASTM D4867](#), Standard Test Method for Effect of Moisture on Asphalt Mixtures