


Image construction for a concave mirror



Physics

Light & Optics

Reflection & refraction of light



Difficulty level

easy



Group size

2



Preparation time

10 minutes




Execution time

10 minutes

This content can also be found online at:

<http://localhost:1337/c/631b118abce9830003710de4>


PHYWE



Teacher information

Application

PHYWE



Reflection at the concave mirror

A concave mirror is a mirror that is curved concavely (inwards).

For example, a concave mirror, which concentrates the light rays in its focal point, is used to harness solar energy.

Pierre Borell, a French physician, first used a concave mirror for medical examination in the 17th century. Slightly modified, they are still used today in medical diagnostics in the form of a forehead mirror,

Other teacher information (1/6)

PHYWE

Prior knowledge



The students should have previously learned the basics of the straight-line propagation of light and the application of the law of reflection. In addition, before carrying out the experiment, it should be explained why objects in our environment can be seen at all and how to construct pixels with the help of main rays.

Principle



Light beams incident on a concave mirror along the optical axis are light beams through the centre of curvature and are therefore reflected into themselves. The "foot points" of the objects "standing" on the optical axis are again imaged on the optical axis. Therefore, only one more image point is sufficient for a clear construction of the image.

Other teacher information (2/6)

PHYWE

Learning objective



With this experiment, the students should learn a procedure that makes it possible to construct the image on the concave mirror with a given object. Selected light beams and their characteristic course are used for this purpose.

Tasks



Why do you see yourself inverted and diminished in a polished spoon?

In this experiment, the students are to investigate this question and examine the formation of images on a concave mirror with the help of selected light beams.

Other teacher information (3/6)

PHYWE

Additional information 1

With this experiment, the student should learn a procedure with the help of which the construction of the image on the concave mirror is possible with a given object. Selected light beams and their characteristic course are used for this purpose.

The experiment is demanding in terms of the students' abilities and experimental skills. However, with careful adjustment and precise experimental work, the experiment is associated with a very high gain in knowledge for the student, especially if corresponding demonstration experiments are carried out with the optical bench as a supplement.

Other teacher information (4/6)

PHYWE

Additional information 2

With this experiment, the essence of the physical experiment can be clearly demonstrated; through the specific specification of the experimental conditions (object width, object size, focal length), one obtains a result with new information content.

By varying the conditions, physical laws can be deduced, but the reverse, deductive path is also possible. The construction of the images with the help of selected light rays and the subsequent experimental verification offer a wide range of possibilities for an interesting lesson design.

The experiment uses an example (geometric constructions) to illustrate the role of mathematics in physical knowledge.

Other teacher information (5/6)

PHYWE

Note

Before carrying out the experiment, explain to the students that objects in our environment are only seen because the light from a light source (e.g. the sun) is reflected by them and enters the eye. It is similar with the formation of an image at the concave mirror. A real image can only ever be created when the light beams emanating from an object point reunite in a point after reflection.

To construct these image points, it is therefore sufficient to select two of the infinite number of light beams and to follow their path from the object point to the image point. For the sake of simplicity, one selects individual light beams ("main beams") and a third light beam serves as a control.

Other teacher information (6/6)

PHYWE

Notes on set-up and procedure

In this experiment, special attention should be paid to the exact position of the light box (e.g. parallel to the optical axis) in the individual experimental steps.

It is therefore advisable to draw a thin auxiliary line beforehand, e.g. parallel to the optical axis. Under these conditions, a high reproducibility of the results can be expected.


Safety instructions

PHYWE

The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Motivation

PHYWE

We encounter mirrors every day in a wide variety of designs. A special type of mirror is the so-called concave mirror. This is a mirror that is curved inwards and concentrates the reflected rays in one point, the focal point.

A typical example of a concave mirror is a spoon, as shown in the picture on the right. Have you ever looked at your reflection in a spoon, do you notice anything special about it?





Image of a spoon on a patterned base

Task

PHYWE



Experimental setup

Why do you see yourself inverted and diminished in a polished spoon?

Investigate the formation of images on a concave mirror with the help of selected light beams.

Equipment

Position	Material	Item No.	Quantity
1	Light box, halogen 12V/20 W	09801-00	1
2	Mirror, concave-convex	09812-00	1
3	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1


Additional equipment

PHYWE

Position	Material	Quantity
1	Ruler (approx. 30cm)	1
2	White paper (DIN A4)	1
3	Circle	1

Set-up

PHYWE



Preparation of the DIN-A4 sheet

Attention!



Make sure that the concave mirror always rests in the middle of the inner curvature on S .

- Prepare a sheet of paper as shown on the left. The distances \overline{FS} and \overline{MS} are 7.2 cm each, the arc of the circle around M has the radius \overline{MS} .
- At a distance of 18 cm from the point S , draw a vertical arrow 2 cm long on the optical axis with a red pencil and label it with G (Subject).
- Draw a thin pencil line exactly parallel to the optical axis as an auxiliary line passing through the tip of the object arrow G .

Procedure (1/6)


PHYWE

- Insert the slit diaphragm into the light box on the lens side. Place the light box and the concave mirror on the sheet of paper.
- Connect the light box to the power supply unit (12 V ~).
- Check the correct position of the concave mirror by first letting the narrow beam of light fall along the optical axis.



Procedure (2/6)

PHYWE




Moving the light box

- Now move the light box until the narrow light beam runs exactly parallel to the optical axis along the auxiliary line and through the tips of the arrow (imaginary object).
- Observe the light beam reflected by the concave mirror and mark the course of the incident and reflected light beam with two crosses. Write down your observations.

Procedure (3/6)

PHYWE



Turning the light box

- Then turn the light box until the light beam passes through the arrowhead of G and the point F (focal point).
- Again observe the light beam reflected by the concave mirror and mark the course of the incident and reflected light beam twice (use a different colour or marker). Write down your observations.

Procedure (3/6)

PHYWE




Picture of the power supply unit

- Switch off the power supply and remove the light box and the concave mirror from the paper.
- Connect the related markings so that the course of the light beams before and after reflection at the concave mirror becomes clear.
- How do the two reflected light beams relate to each other? Write down your observations.

PHYWE

Report



Task 1

PHYWE

Consider the following statement:

Both light beams intersect after reflection at a point below the optical axis.

☐ True☐ False☒ Check

Task 2

PHYWE

Why do you see yourself inverted and diminished in a polished spoon?

A polished spoon is a small . The and the laws governing the of selected light beams thus apply to it. Since the is outside the focal length, the result is a , inverted image.

observer

concave mirror

double

reduced

course

law of reflection

☒ Check

Task 2

PHYWE

Why do you see yourself inverted and diminished in a polished spoon?

A polished spoon is a small . The and the laws governing the of selected light beams thus apply to it. Since the is outside the focal length, the result is a , inverted image.

observer

concave mirror

double

reduced

course

law of reflection

☒ Check

Task 3

PHYWE

Draw a straight line from the tip of the object arrow G through the point M to the mirror (centre beam).

What do you see?

The also passes through the of the other two selected rays and thus through the of the image arrow. A light beam that falls through the onto the concave mirror is reflected (angle of incidence = = 0°). If it touches the tip of the object on its way or is emitted by it, it must also touch of the image, since light beams emitted from a point on the object contribute to the formation of the image.

centre ray

the tip

centre of curvature M

intersection

in itself

all

angle of reflection

tin

Task 4

PHYWE

Why is it sufficient to draw two excellent rays of light from objects "standing" on the optical axis, only from the tip of the object, in order to construct the image?

Light beams incident of the optical axis on a concave mirror are light beams through the centre of curvature and are therefore reflected parallel displaced. The "foot points" of the objects "standing" on the optical axis become the arrowheads and vice versa. Therefore, only one more image point is sufficient for a clear construction of the image.

Light beams incident on a concave mirror along the optical axis are light beams through the centre of curvature and are therefore reflected into themselves. The "foot points" of the objects "standing" on the optical axis are again imaged on the optical axis. Therefore, only one more image point is sufficient for a clear construction of the image.

Additional task

PHYWE

The lens grinder equation applies to the image formation at the concave mirror:

$$\frac{1}{f} = \frac{b}{g} + \frac{1}{g} ,$$

where $f = 7,2 \text{ cm}$ = focal length, b = image width = distance BS, g = object width = distance GS .

Check your measurement results in table 1 with this equation.