

# Illuminance (inverse square law)



Physics

Light &amp; Optics

Dispersion of light



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:



<http://localhost:1337/c/5f4e975b38db8d0003265b8e>

PHYWE

## Teacher information



## Application

PHYWE



range of an incandescent lamp

How far does the light of an incandescent lamp actually reach? A flashlight? Or a headlight?

This experiment is about the decrease of illuminance with increasing distance from the light source.

## Other teacher information (1/3)

PHYWE

## Prior knowledge



The pupils should first know the basics of linear propagation of light and be able to use a light box, and in connection with the measurement of the side lengths of the illuminated surfaces, the pupils must draw on their knowledge of the core and penumbra. Because the light source is not punctiform, the edges of the shadows are not sharp.

## Scientific principle



The area illuminated by a light box is determined depending on the distance to the light box.

## Other teacher information (2/3)

PHYWE

## Learning objective



This experiment is suitable both as an introductory experiment and as a confirmation experiment for the validity of the law of distance.

## Tasks



Investigations such as illuminance  $E$  with the distance  $r$  between a light source and an illuminated surface.

## Other teacher information (3/3)

### Notes on construction and implementation

The teacher should make sure that the graph paper on the screen is not wavy but tightly stretched and that when marking the edges of the illuminated areas, the screen must not be shifted or tilted, so that the measuring errors do not become too large

## Safety instructions



- Halogen lamps become warm during prolonged use
- Avoid looking directly into the light source



## Safety instructions


PHYWE



- Halogen lamps become warm during prolonged use
- Avoid looking directly into the light source

PHYWE

## Student Information



## Motivation

PHYWE



### On the road at night

How far can you see at night on the road when you are riding your bike?

From what distance can others recognize you when you have your lamp on?

The answers to these questions are hidden in the distance dependence of illuminance.

## Equipment

Position	Material	Item No.	Quantity
1	Optical profile-bench for student experiments, l = 600 mm	08376-00	1
2	Light box, halogen 12V/20 W	09801-00	1
3	Bottom with stem for light box	09802-20	1
4	Diaphragm with square	09816-03	1
5	Slide mount for optical bench	09822-00	1
6	Screen, white, 150x150 mm	09826-00	1
7	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1

## Equipment

PHYWE

Position	Material	Item No.	Quantity
1	<a href="#">Optical profile-bench for student experiments, l = 600 mm</a>	08376-00	1
2	<a href="#">Light box, halogen 12V/20 W</a>	09801-00	1
3	<a href="#">Bottom with stem for light box</a>	09802-20	1
4	<a href="#">Diaphragm with square</a>	09816-03	1
5	<a href="#">Slide mount for optical bench</a>	09822-00	1
6	<a href="#">Screen, white, 150x150 mm</a>	09826-00	1
7	<a href="#">PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A</a>	13506-93	1

## Set-up (1/5)

PHYWE



Assembly of the optical bench

Build up the optical bench from the two tripod rods and the variable tripod foot and place the scale on the front tripod rod.



## Set-up (2/5)

PHYWE



Using the light box with base and handle

- Place the base with handle under the light box.

## Set-up (3/5)

PHYWE




Positioning of the light box

- Clamp it into the left part of the tripod base with the lens side facing away from the optical bench.
- Slide an opaque shade in front of the lens and the shade with the square opening into the shaft at the other end of the lamp.

## Set-up (4/5)

PHYWE




Attaching the millimetre paper

- Place the graph paper sheet on the screen, fold the protruding parts backwards and fasten the paper with 3 paper clips tightly to the screen.

## Set-up (5/5)

PHYWE



Attaching the millimetre paper

- Place the shade on the optical bench using the tab by the lamp.

## Procedure (1/2)

PHYWE




Connecting the light box

- Connect the lamp to the power supply unit (12 V~) and switch it on.

## Procedure (2/2)

PHYWE



Positioning of the shield

- Slowly move the screen to the right while observing the illuminated area.

## Procedure (2/2)

PHYWE




Marking of the illuminated area

- Place the screen at the 6 cm mark. (The screen has now a distance of  $r = 6$  cm from the light source (filament of the lamp in the lamp)).
- Mark the edges of the illuminated area with a pencil (dashed).

## Procedure (2/2)

PHYWE




Marking of the illuminated area

- Place the screen at 12 cm, 18 cm and 24 cm intervals.
- Mark the respective edges of the illuminated area in the same way.
- Switch off the power supply unit.
- Loosen the graph paper from the screen and draw the illuminated areas by joining the marks of their edges.
- Measure the side lengths  $a$  and  $b$  of the illuminated areas and enter them in the table in the protocol.

PHYWE

# Report



## Task 1

PHYWE

Change of the illuminated area with increasing distance of the screen from the light source: Enter your measured values in the table. For the illuminated areas, refract the areas  $A = a \cdot b$  and the squares of their respective distances  $r$  from the light source.

r in cm	a in cm	b in cm	a*b in cm <sup>2</sup>	r <sup>2</sup> in cm <sup>2</sup>	A/r <sup>2</sup>
3 (aperture)					
6					
12					
18					
24					

## Task 2

PHYWE

Compare the values of  $A = a \cdot b$  and  $r^2$ . What is probably the relationship between  $A$  and  $r^2$ ?

$A$  and  $r^2$  are (probably)  to each other.

☒ Check

## Task 3

PHYWE

Check your guess: Calculate the quotients  $A/r^2$  (to 2 decimal places) and enter your results in the last column of Table 1 in the section "Observations and Measurement Results". Write down the relationship between  $A$  and  $r^2$  in mathematical form and in words.

$A / r^2 =$

.

$A \sim r^2$ , i.e.  of the illuminated area and  of its distance from the light source are proportional to each other.

☒ Check

## Task 4

PHYWE

The light emitted by the light source is therefore distributed over four times the area of the surface, e.g. at double the distance  $r$ , i.e. illuminance  $E$  (or brightness) has fallen to a quarter.

So what is the relationship between  $E$  and  $r$ ?

## Task 5

PHYWE

The illuminance  $E$  is also proportional to the luminous intensity  $I$  of the light source.

What is the relationship between  $E$ ,  $I$  and  $r$ ?