

# Light path and focal length of a concave lens



Physics

Light &amp; Optics

Reflection &amp; refraction of light



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/5f4e5c4938db8d0003265980>

PHYWE



## Teacher information

### Application

PHYWE



single-lens imaging

We all use optical devices every day. These are mobile phone cameras, cameras, microscopes, binoculars and many, many more.

They all use lenses to enable the imaging of objects. Often these are combinations of convex and concave lenses.

This experiment deals with the imaging properties of concave lenses and thus strengthens the understanding of optical devices.

## Other teacher information (1/4)

PHYWE

### Prior knowledge



Students should have previously learned the basics of linear propagation of light and the concepts of diffraction and refraction.

### Scientific principle



In connection with the experiments on the refraction of light, this experiment is of particular importance. The knowledge of the law of refraction is consolidated and transferred to a new subject. The experiments on the refraction of light on concave lenses serve to consolidate the knowledge and experimental skills gained in the experiments on convex lenses.

## Other teacher information (2/4)

PHYWE

### Learning objective



The aim of this experiment is to observe the course of light as it passes through concave lenses and to determine the focal point and thus the focal length using the known methods. On the other hand, however, the position of the (virtual) focal point in front of the lens reveals a significant difference to convex lenses and thus prepares the introduction of the term "virtual image".

### Tasks



Examination of the course of light through a plano-concave lens and determination of the focal length.

## Other teacher information (3/4)

The term "diverging lenses" is commonly used for lenses that are thicker at the edges than in the centre. However, the condition under which this statement is correct is very often neglected. Concave air-filled lenses in water - i.e. diverging lenses starting from the outer shape - have a collecting effect.

Correspondingly, convex air-filled lenses in water - i.e. convergent lenses starting from the outer shape - have a diverging effect.



## Other teacher information (4/4)

PHYWE

### Instructions for construction and implementation

Care must be taken to ensure that the adjustment of the concave lens (flat surface at the vertical of the line cross, unbroken course of a light beam incident along the optical axis) is carried out very carefully by the student in order to achieve a clear and convincing experimental result.

If necessary, difficulties for the students could arise from the fact that a virtual focal point occurs in the concave lens as in the domed mirror. This means that a rearward extension of the refracted, diverging light bundles must be made.

## Safety instructions


PHYWE



- Halogen lamps become warm during prolonged use
- Avoid looking directly into the light source

PHYWE

## Student Information



## Motivation

PHYWE



Binoculars


### Optical devices:

We use binoculars, mobile phone cameras or microscopes frequently, if not daily, without thinking about what is actually inside.

In optical devices these are usually combinations of different lenses with different optical properties. One of these lenses and its imaging properties will be explained in this experiment.

## Tasks

PHYWE



Experiment set-up

### What distinguishes concave lenses from convex lenses?

1. Examine the course of light through a plano-concave lens and determine the focal length.

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Light box, halogen 12V/20 W</a>	09801-00	1
2	<a href="#">Block, planoconcave lens, fl-100mm</a>	09810-05	1
3	<a href="#">PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A</a>	13506-93	1


## Additional equipment

PHYWE

Position	Material	Quantity
1	White paper (DIN A4)	1
2	Ruler (approx. 30 cm)	1

## Set-up (1/2)

PHYWE



Production of the line cross

### Look out!

Make sure that in all experiments with the flat surface the lens is exactly on the vertical line of the line cross and that the model body does not change its position when moving the light box.


### 1. course of the light through a plano-convex lens.

Draw a right-angled line cross in the middle of a sheet of paper. The point of intersection of the lines is *M*. Draw in each 3 cm distance from *M* one mark on each vertical line.



## Set-up (2/2)

PHYWE




Setting up the light box

- Connect the light box to the power supply unit (12 V ~)
- Place the light box with the lens side, but without the aperture, on the edge of the sheet.
- Place the plano concave lens (roughened side down) with the flat surface exactly on the vertical line of the line cross within the two marks.

## Procedure (1/5)

PHYWE




Illumination of a concave lens

- Observe the course of the parallel light as it passes through the lens.

## Procedure (2/5)

PHYWE




Lateral incidence of light into a concave lens

- Now move the light box to the positions according to the illustrations and observe the course of the light again.
- Note your observations in the protocol.

## Implementation (3/5)

PHYWE




Using the light box with five slit diaphragm

- Insert the five slit diaphragm into the light box on the lens side and place it at a distance of about 10 cm from the inwardly curved (concave) surface of the model body.
- The central light beam should be incident exactly along the optical axis. If it is not along the optical axis, carefully move the lens slightly along the vertical line.

## Procedure (4/5)

PHYWE




Marking the position of the lens

- Mark the outlines of the lens with a thin pencil stroke.
- Describe the course of the narrow light beams as they pass through the lens, in particular the light path within the lens in the protocol.

## Procedure (5/5)

PHYWE



Sketch of the beam path

- Always mark with two crosses the course of the light beams above and below the optical axis before and after passing through the lens.
- Switch off the power supply and remove the light box and model body from the paper.
- Connect the matching crosses so that the course of the light beams before and after passing through the lens and, after appropriate connection, also within the lens is visible.

PHYWE



# Report

## Observation

PHYWE

Write down your observations.

a) Observation of the course of light without aperture:

The light incident in parallel on a concave lens is  by the lens and diverges behind the lens, it diverges.

b) Observation of the course of light with a five slit diaphragm:

The narrow beams of light incident on a concave lens parallel  are  as they pass through the glass of the lens and diverge behind the lens.

☒ Check

## Task 1

PHYWE

Using your observations, describe how parallel light behaves when incident on a plano-concave lens.

Parallel to the ,  incident on a   
run outwards after the , they .

diverge

concave lens

light beams

fraction

optical axis

☒ Check

## Task 2

PHYWE

Extend the refracted light rays with a coloured pencil until they intersect. What can you determine about the position of the intersection?

The  of the rear extensions of the   
is in front of the  (on the side of the light incidence) on the  
.

concave lens

intersection

optical axis

refracted light beams

☒ Check

## Task 3

PHYWE

Designate the point of intersection of the refracted rays with F'.

Measure the distance of the point F' (the focal point) from the center M and enter the value here.

$\overline{MF'}$  =  cm

Note: Since the lens can no longer be considered a thin lens, the point of intersection of the outer rays is slightly closer to the lens than that of the inner rays. For the distance  $\overline{MF'}$  the average value should be given.

☒ Check

## Task 4

PHYWE

Why are light bundles incident on a concave lens along the optical axis not refracted?

For a light beam incident on a concave lens along the , the  and thus the  is 0°, there is .

☒ Check

## Task 5

PHYWE

What distinguishes concave lenses from convex lenses?

In contrast to convex lenses, parallel incident light is not [ ] at a concave lens, but [ ]. The ( [ ]) focal point of concave lenses is [ ] of the lens, whereas with convex lenses the ( [ ]) focal point is [ ] the lens (viewed from the side of the light incidence). Lenses that have this scattering property are also called [ ].

collected

virtual

diffusing lenses

real

scattered

in front

behind

 Check

Slide	Score / Total
Slide 21: Light pattern	0/3
Slide 22: Parallel light on plano-concave lens	0/5
Slide 23: Position of the intersection point	0/4
Slide 24: Point of intersection of the refracted rays	0/1
Slide 25: Light beam along the optical axis	0/4
Slide 26: Difference between concave and convex lenses	0/7

Total amount  0 / 24 Solutions Repeat