


Law of imagery for a convex lens



Physics

Light & Optics

Optical devices & lenses



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/62dbc704a52f910003dffb87>


PHYWE



Teacher information

Application

PHYWE



Experimental setup

Convex lenses, also called converging lenses, can produce a magnified image. They are an important element of ray optics and are therefore widely used in optical instruments and in photographic lenses.

Other teacher information (1/4)

PHYWE

Principle



Incident light that is parallel to the optical axis is focused by the convex lens at the focal point. This can produce a magnified real image.

Learning objective



The students are asked to observe the optical effect of a convex lens and to explain the imaging law $1/f = 1/b + 1/g$

Other teacher information (2/4)

PHYWE

Task



- The students should investigate what relationship there is between the focal length f , the subject matter g and the image width b when real images are created with the help of a convex lens.
- For this purpose, for given object widths g the image width b will be measured and noted in a table in the report.

Other teacher information (3/4)



- The law of mapping is generally first derived theoretically or simply given.
- It would also be conceivable for the students to experiment without prior knowledge of the law of mapping. In this case, however, the student would not understand why the reciprocals of f , g and b are to be formed.

Other teacher information (4/4)

PHYWE

Notes on set-up and procedure

- The proposed experimental steps are intended to ensure that both enlarged and reduced images of the same size are produced.
- As another experiment technique, the teacher can also suggest that the students first place the screen in different places and then move the lens until the image is sharp.
- Do not carry out the experiment with the convex lens $f = 50 \text{ mm}$ because the enlarged images created with it can no longer be focused. The condition that the rays are close to the axis is then no longer fulfilled.

Safety instructions


PHYWE



- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Motivation

PHYWE




Magnifying glass as an example of a convex lens

Convex lenses, also called converging lenses, can produce a magnified image. They are a central component of optics and are therefore often found in everyday devices, such as telescopes, camera lenses or even eyeglasses.

How do convex lenses work?

Tasks

PHYWE



Experimental setup

- Investigate the relationship between the focal length f , the subject matter g and the image width b that exists when real images are created with the help of a convex lens.
- For this purpose, for the given object widths g you shall measure the image width b and note in a table in the report.

Equipment

Position	Material	Item No.	Quantity
1	Optical profile-bench for student experiments, l = 600 mm	08376-00	1
2	Light box, halogen 12V/20 W	09801-00	1
3	Bottom with stem for light box	09802-20	1
4	Lens on slide mount, f=+100mm	09820-02	1
5	Slide mount for optical bench	09822-00	2
6	Screen, white, 150x150 mm	09826-00	1
7	Object -L-, glass bead	11609-00	1
8	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1

Set-up (1/3)

PHYWE

- Assemble the optical bench from the two tripod rods and the variable tripod foot and place the scale on the front tripod rod.
- Place the base with stem under the light box.



Set-up (2/3)


PHYWE

- Clamp the light box in the left part of the tripod base so that the lens side faces away from the optical bench.
- Slide an opaque shade in front of the lens and the Perl-L into the shaft at the other end of the luminaire.



Set-up (3/3)

PHYWE




Experimental setup

- Place the lens and, with the help of the slide, the screen on the optical bench.

Procedure (1/3)


PHYWE




- Connect the lamp to the power supply unit (12 V~) and switch it on.
- Place the lens at a distance of about 20 cm from the Perl-L and move the screen until the image of the Perl-L is as sharp as possible.

Procedure (2/3)

PHYWE




- Measure the distance of the lens from the screen, the image width b and carry for the object width $g = 200 \text{ mm}$ in Table 1 in the report.
- Move the lens to the left and determine for two different object widths ($g > 150 \text{ mm}$) the image widths.



Procedure (3/3)


PHYWE



Displacement of the lens

- Then move the lens to the right and also determine for two different object widths ($g > 200 \text{ mm}$) the image widths.
- Enter the values for g and b in Table 1 in the report.
- Switch off the power supply unit.

PHYWE



Report

Table 1

PHYWE

Enter your measurements in the table.

g in mm b in mm $1/g$ in mm $1/b$ in mm $1/f$ in mm $1/g + 1/b$ in mm

200	0.005	0.01			

Task 1

PHYWE

Compare the values in the last two columns of the table. What do you find?

- ☐ The values contained in the last two columns are the same except for minor deviations.
- ☐ The values contained in the last column are all greater than the values from the penultimate column.
- ☐ No correlation can be found between values in the last and penultimate columns.

☒ Check

Task 2

PHYWE

Calculate the mean value for the sum $1/g + 1/b$ and compare this with the value $1/f$. What is the result in mathematical form?

- ☐ $1/f < 1/g + 1/b$
- ☐ $1/f = 1/g + 1/b$
- ☐ $1/f = 1/g = 1/b$

☒ Check

Task 3

PHYWE

Why was it recommended in the experimental description to choose $g > 150 \text{ mm}$?

- ☐ If $g < 150 \text{ mm}$ the image on the screen is only a virtual image.
- ☐ If $g < 150 \text{ mm}$ the image is created at such a great distance that the optical bench is no longer sufficient to catch it on the screen.
- ☐ If $g < 150 \text{ mm}$ the image on the screen is merely reversed.

☒ Check

Task 4

PHYWE

What would happen in the event of $g = 100 \text{ mm}$?

- ☐ $b = 0$ i.e. the image width is equal to zero.
- ☐ $b = g$ i.e. the object width is equal to the image width.
- ☐ $1/b = 0$ i.e. the reciprocal of the image width is equal to zero.
- ☐ $1/f = 1/g$ i.e. the object distance is as large as the focal length.
- ☐ $b = \infty$ i.e., the image lies in the infinite.

☒ Check

Task 5

PHYWE

Fill in the blanks: What would happen in the event of $g < 100 \text{ mm}$?

For $g < 100 \text{ mm}$ the object distance is than the focal length f .  images are created. These are not catchable on the .

☒ Check


Magnifying glass as an example of a convex lens

Slide	Score/Total
Slide 20: Comparison of the values	0/1
Slide 21: mathematical relation	0/1
Slide 22: Case $\backslash g > 150 \backslash, \backslash \mathrm{mm} \backslash$	0/1
Slide 23: Case $\backslash g = 100 \backslash, \backslash \mathrm{mm} \backslash$	0/3
Slide 24: Case $\backslash g < 100 \backslash, \backslash \mathrm{mm} \backslash$	0/3

Total  ★ 0/9

 Solutions

 Repeat

 Export text