

What does the LED spectrum of the transmission grating look like?



Physics

Light & Optics

Spectrometry & refractometry



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/5f5060de37ffe20003f0fff9>

PHYWE

Teacher information



Application

PHYWE



Experiment set-up

What does the LED spectrum of the transmission grating look like?

The wavelength of light can be determined in many different ways.

In this experiment the students learn the so-called objective method with a transmission grating.

The name "transmission grating" means that the light passes through the grating and interferes with it.

Other teacher information

PHYWE

Notes



Since scattered light hardly affects this measurement, the classroom only needs to be slightly darkened so that the interference is visible on the screen. Distances and distances must be measured very accurately, since even small inaccuracies can cause large deviations in the result.

Task



Determine the wavelength of maximum intensity with a transmission grating.

Safety instructions


PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student Information



Motivation

PHYWE




Experiment set-up

What does the LED spectrum of the transmission grating look like?

The wavelength of light can be determined in many different ways.

In this experiment you will learn the so-called objective method with a transmission grating.

The name "transmission grating" means that the light passes through the grating and interferes with it.



Equipment

Position	Material	Item No.	Quantity
1	Support base, variable	02001-00	1
2	Support rod, stainless steel, l = 600 mm, d = 10 mm	02037-00	2
3	Slide mount without angle scale	09851-02	1
4	Diaphragm holder, attachable	11604-09	2
5	Lens on slide mount, f=+100mm	09820-02	1
6	Screen, semitransparent, 150x150mm ²	09851-03	1
7	Lens on slide mount, f=+300mm	09820-04	1
8	Grating, 500 lines/mm, in slide frame, glassless	09851-16	1
9	LED - red, with series resistor and 4 mm plugs	09852-20	1
10	Stray light tube for LED, Di = 8 mm, l = 40 mm	09852-01	1
11	Measuring tape, l = 2 m	09936-00	1
12	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
13	Connecting cord, 32 A, 750 mm, red	07362-01	1
14	Connecting cord, 32 A, 750 mm, blue	07362-04	1

Set-up (1/2)

PHYWE

- Connect the stray light tube to the LED according to the illustrations.
- Place the components on the stand material.
- The illumination gap and the grid are not yet needed.



Set-up (2/2)


PHYWE

- The LED is connected to the power supply unit (make sure the polarity is correct).
- The power supply unit is set to 6 V.



Procedure (1/2)



PHYWE



- The lens is moved back and forth on the tripod material.
- As long as a sharp (and small) light spot can be seen up to the screen.
- The illumination gap and the grating are placed together in an aperture holder on the rider with the lens facing the screen.

Procedure (2/2)

PHYWE



- An interference image is now visible on the screen.
- Measure the distance between the first two maxima.
- Note this as $2 \cdot a$
- Measure the distance between the grid and the screen.
- Note this as l .

PHYWE



Report

Task 1

PHYWE

Note down the measured values

for a and l .

$a =$ cm

The grating has 500 lines per mm and the wavelength is 632nm

Combine both equations (right) and solve them according to λ at

$\lambda =$ \cdot $\arctan(\frac{\text{}}{\text{}})$

l $\sin($ a g

☒ Check

What value do you get for λ ?

$\lambda =$

The following equations/formulas are known:

$$\sin(\alpha) = \frac{\lambda}{g}$$

$$\tan(\alpha) = \frac{a}{l}$$

λ = wavelength of light

g = the grid constant

a = distance to the first interference maximum

l = distance between grille and screen

Task 2

PHYWE



Combine both equations (right) and solve them according to λ at

By solving the formula $\tan(\alpha) = \frac{a}{l}$ to a and inserting it into the formula $\sin(\alpha) = \frac{\lambda}{g}$ which are looking for l is dissolved, the formula is obtained: $\lambda =$

. By inserting the values of a , l and g (500 strokes per mm) you get $\lambda =$ \times

$=$. The wavelength of the red LED is specified as 632 nm, which means that there is a deviation of just under in this measurement.

Slide

Score/Total

Slide 13: Formula

0/4

Slide 14: Solving the formula

0/5

Total amount