


Generation of an induced voltage with permanent magnets



Physics

Electricity & Magnetism

Electromagnetism & Induction

Physics

Electricity & Magnetism

Electric generator, motor, transformer



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:

<http://localhost:1337/c/5f4fde2137ffe20003f0fe6b>


PHYWE

Teacher information



Application

PHYWE



Experiment set-up

Electromagnetic induction, or Faraday induction, has many possible applications and it is impossible to imagine our everyday life without it.

For example, microphones and loudspeakers convert acoustic signals into electrical signals and vice versa.

Generators such as the bicycle dynamo, water turbines on dams or wind turbines also use the principle of induction to produce electricity.

Other teacher information (1/3)

PHYWE

Prior knowledge



Students should know that a conductor carrying current is surrounded by a magnetic field. They should also know that this magnetic field can cause mechanical movement in interaction with another magnetic field. Here they should demonstrate the reversal of this phenomenon.

Scientific principle



The change of the magnetic field on a conductor (here: a coil) causes an electric field. The law of induction states

$$\text{rot} \vec{E} = \vec{\nabla} \times \vec{E} = -\frac{\partial \vec{B}}{\partial t}$$

Other teacher information (2/3)

PHYWE

Learning objective



The students should learn that electrical energy can be generated by mechanical movement. They should recognize that a voltage is only induced when the magnetic field within a coil (or conductor loop) changes.

Tasks



The students are to set up a galvanometer to measure induced voltages and generate an induction voltage by the interaction of a permanent magnet with a coil.

Other teacher information (3/3)

PHYWE

Notes

The galvanometer has the advantage over the multiple measuring instrument that it can deflect to both sides and therefore allows the indirect observation of currents of different directions.

Quantitative statements on the induction voltage are not aimed at with this experiment.

As the measuring mechanism of the galvanometer also has a permanent magnet, care should be taken that the bar magnet to be experimented with is not moved in the immediate vicinity of the galvanometer. Otherwise the sensitive measuring mechanism will react to the magnet and falsify the measured values to be expected by induction.


Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE



Student Information

Motivation

PHYWE



Wind farm with wind turbines

If one wants to convert mechanical energy into electrical energy, this is done primarily through electromagnetic induction.

You probably know this principle from applications such as the bicycle dynamo for the Farad lamps or from wind turbines in wind parks or water turbines in dams. In all these applications, magnets are usually moved relative to conductor loops, which induces an electrical voltage.

In this experiment you will investigate how an electric voltage can be induced with the help of a permanent magnet.

Tasks

PHYWE



- Set up a galvanometer for qualitative voltage measurement.
- Move a permanent magnet in various ways relative to a coil connected to the galvanometer.

Equipment

Position	Material	Item No.	Quantity
1	Straight connector module with socket, SB	05601-11	1
2	Coil, 400 turns	07829-01	2
3	Iron core, I-shaped, laminated	07833-00	1
4	magnet, l = 72mm, rodshaped, colored poles	07823-00	1
5	Galvanometer movement	07875-00	1
6	Galvanometer scale	07876-00	1
7	Notch bearing with plug	07877-00	1
8	Connecting cord, 32 A, 500 mm, red	07361-01	2

Set-up (1/2)

PHYWE


Set up the galvanometer according to the illustrations. Assemble the coil with 400 turns, the chine bearing, the galvanometer scale and finally the galvanometer movement.

Plug the finished galvanometer into a cable module with socket.



Set-up (2/2)

PHYWE



Complete the test setup as shown in the adjacent figure. Connect the second coil with the two connecting leads to the coil of the galvanometer.

Position the second coil as far away as possible from the galvanometer so that only the induced voltage affects the galvanometer and not the magnetic interaction of the permanent magnet with the galvanometer movement.

Procedure

PHYWE

Carry out the following test steps one after the other. Observe the deflection of the galvanometer pointer and write down your observations.

1. Move the magnet with the **North Pole** into the coil and out again after a short pause.
2. Move the magnet with the **South Pole** into the coil and out again after a short pause.
3. Move the magnet in and out of the coil faster.
4. Move the coil towards the magnet and after a short pause move it away from the magnet again.
5. Let the magnet rest in the coil.
6. Rotate the magnet by tilting it slightly around its longitudinal axis without moving it.

Note The movements in steps 1, 2 and 4 should be as fast as possible.

PHYWE

Report




Table 1

PHYWE

Enter your observations for the pointer deflections (to the left/right or smaller/bigger) in the table.

Movement	1a. North Pole into Coil	1b. North Pole out of Coil	2a. South pole into coil	2b. South pole out of coil	3. faster movement of the magnet
deflection of the pointer					

Movement	4a. Coil towards the	4b. Coil away from magnet	5. magnet rests in the coil	6. rotation of the magnet around the longitudinal axis
deflection of the pointer				

Task 1

PHYWE

The voltage that the galvanometer indicated during the test is called induction voltage. The process by which it is generated is called electromagnetic induction.

What is the direction of the induction voltage dependent on (compare the results of steps 1. and 2.)?

- ☐ The direction depends on which pole of the magnet faces the coil.
- ☐ The direction depends on whether the movement of the magnet is into or out of the coil.
- ☐ The direction depends on how fast the magnet is moved.

✓ Check

Task 2

PHYWE

What does the strength of the induction voltage depend on (compare the result of step 3. with the previous ones)?

- ☐ The strength of the induction voltage depends on the speed of the movement.
- ☐ The strength of the induction voltage depends on the pole of the magnet.
- ☐ The strength of the induction voltage depends on the direction of movement of the magnet.

☒ Check

Task 3

PHYWE

Complete the following text. Drag the words to the appropriate places.

For the generation of the it whether the moves e.g. towards the coil or the towards the magnet.

Not required:

coil

doesn't matter

important

induction voltage

magnet

☒ Check

Task 4

PHYWE

Complete the following definition for induction voltages. Drag the words to the appropriate places.

A is when the enclosed by magnetic field the changes.

magnetic field

voltage

induced

induction coil

 Check

Slide

Score/Total

Slide 16: Comparison of the measurements - direction of the induced...

0/2

Slide 17: Comparison of the measurements - strength of the inductio...

0/1

Slide 18: Result of the relative movement

0/5

Slide 19: Definition of induction

0/4

Total amount

  0/12

Solutions



Repeat



Exporting text