


At what wavelength does an LED light up?



Physics

Modern Physics

Solid state physics



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/615c3b5dbd5b06000320c685>


PHYWE

Teacher information



Application

PHYWE



Test setup

The wavelength of light influences the results of many experiments. Thus, by knowing the wavelength of emitted light, much can be said about the nature of different materials.

Determining the wavelength of light is a good way to check the understanding of diffraction at the grating, as this is necessary to determine the wavelength.

Other teacher information (1/2)

PHYWE

Previous



Students should be familiar with diffraction at the grating.

Principle



Light rays hit the grating and are diffracted there with different degrees of rigidity according to their wavelength. Thus, the position of the first secondary maximum with a known grating constant depends only on the distance of the grating from the screen and the wavelength of the light used.

Other teacher information (2/2)

PHYWE

Learning



The aim of this experiment is to consolidate the students' understanding of diffraction at the grating.

Tasks



- Measurement of the first secondary maximum and determination of the wavelength of the light used.

PHYWE




Student Information

Motivation

PHYWE

The determination of the wavelength of light is a vast field, which is especially applied in astronomy and materials research.

Whether it is used in spectroscopy to find out what the atmosphere of a distant planet is made of, or to find out how fast a car fell into a speed trap, the determination of the wavelength of the emitted light accompanies many areas of life.



Equipment

| Position | Material | Item No. | Quantity |
|----------|---|----------|----------|
| 1 | LED - blue, with series resistor and 4 mm plugs | 09852-40 | 1 |
| 2 | LED - green, with series resistor and 4 mm plugs | 09852-30 | 1 |
| 3 | LED - red, with series resistor and 4 mm plugs | 09852-20 | 1 |
| 4 | LED - UV, with series resistor and 4 mm plugs | 09852-50 | 1 |
| 5 | Grating, 500 lines/mm, in slide frame, glassless | 09851-16 | 1 |
| 6 | Support base, variable | 02001-00 | 1 |
| 7 | Support rod, stainless steel, l = 600 mm, d = 10 mm | 02037-00 | 2 |
| 8 | Slide mount without angle scale | 09851-02 | 2 |
| 9 | Diaphragm holder, attachable | 11604-09 | 2 |
| 10 | Measuring tape, l = 2 m | 09936-00 | 1 |
| 11 | PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A | 13506-93 | 1 |
| 12 | Connecting cord, 32 A, 750 mm, red | 07362-01 | 1 |
| 13 | Connecting cord, 32 A, 750 mm, blue | 07362-04 | 1 |

Structure (1/2)

PHYWE

- Note: To perform this experiment, it is advantageous if the room is only slightly darkened and the measurements are not performed directly against daylight.
- Set up the experiment according to steps 1 to 6 in the order LED, measuring scale, grid on the stand material.



Step 1



Step 2




Step 3

Structure (2/2)

PHYWE


- The prepared measuring scale is plugged onto the diode so that the head of the LED just sticks out. Pay attention to the correct polarity when connecting the LED!



Step 4



Step 5




Step 6

Procedure (1/2)

PHYWE

- Select the LEDs with the colors UV, blue, red and green one after the other.
- After connecting an LED, set a voltage of 6 V so that the LED lights up clearly.
- The eye is used to get very close to the grid from behind.
- Move the grid forward/backward until the first maximum is on the middle of a thick line.





The accuracy is increased by setting a distance of 8 or 10 cm for the first maximum.


Procedure (2/2)

PHYWE

- Note the distance of the first maximum to the center as a_k in Table 1 in the Protocol.
- Measure the distance from the grid to the tip of the LED and note it as distance l in Table 1.
- Repeat the experiment for the other diodes.



PHYWE



Report

Task 1

PHYWE

| LED color | a_k in cm | l in cm | λ in nm |
|-----------|----------------------|----------------------|----------------------|
| red | <input type="text"/> | <input type="text"/> | <input type="text"/> |
| green | <input type="text"/> | <input type="text"/> | <input type="text"/> |
| blue | <input type="text"/> | <input type="text"/> | <input type="text"/> |
| UV | <input type="text"/> | <input type="text"/> | <input type="text"/> |

Task 2

PHYWE

What equations are needed to determine the lattice constant?

☐ $\sin(a_k) = k \cdot \lambda / g$

☐ $\tan(a_k) = a_k N$

☐ $\cos(a_k) = k \cdot \lambda / g$

☒ Check

Task 3

PHYWE

What formula does this give for the lattice constant?

☐ $g = k \cdot \sin(a_k N)$

☐ $g = k \cdot (\lambda / \cos(\arctan(a_k N)))$

☐ $g = k \cdot (\lambda / \sin(\arctan(a_k N)))$

| Slide | Score / Total |
|---|---------------|
| Slide 14: Determination of the lattice constant | 0/2 |
| Slide 15: Equation 2 | 0/1 |

Total score  0/3



Show solutions



Repeat



Export text