

Band gap of germanium



Physics

Modern Physics

Solid state physics



Difficulty level

hard



Group size

2



Preparation time

45+ minutes



Execution time

45+ minutes

This content can also be found online at:

<http://localhost:1337/c/600a9f516080870003ab8bb2>

PHYWE



General information

Application

PHYWE



Fig.1: Experimental set-up for the determination of the band gap of germanium

Semiconductors have wide applications in electric circuits. Knowledge about the band gap of the used semiconductor material is very important to use them effectively.

Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

The conductivity of a germanium testpiece is measured as a function of temperature. The energy gap is determined from the measured values.

Other information (2/2)

PHYWE

**Learning
objective****Tasks**

The goal of this experiment is to determine the band gap of germanium.

1. The current and voltage are to be measured across a germanium test-piece as a function of temperature.
2. From the measurements, the conductivity σ is to be calculated and plotted against the reciprocal of the temperature T . A linear plot is obtained, from whose slope the energy gap of germanium can be determined.

Theory (1/2)

PHYWE

The conductivity σ is defined as following:

$$\sigma = \frac{1}{\rho} = \frac{l \cdot I}{A \cdot U} \left[\frac{1}{\Omega \text{m}} \right]$$

with ρ = resistivity, l = length of test specimen, A = cross section, I = current, U = voltage. (Dimensions of Ge-plate $20 \times 10 \times 1 \text{ mm}^3$)

The conductivity of semiconductors is characteristically a function of temperature. Three ranges can be distinguished: at low temperatures we have extrinsic conduction (range I), i.e. as the temperature rises charge carriers are activated from the impurities. At moderate temperatures (range II we talk of impurity depletion, since a further temperature rise no longer produces activation of impurities. At high temperatures (range III), it is intrinsic conduction which finally predominates (see Fig. 2). In this instance charge carriers are additionally transferred by thermal excitation from the valence band to the conduction band. The temperature dependence is in this case essentially described by an exponential function.

Theory (2/2)

PHYWE

$$\sigma = \sigma_0 \cdot \exp\left(-\frac{E_g}{2kT}\right)$$

(E_g = energy gap, k = Boltzmann's constant, T = absolute temperature).

The logarithm of this equation

$$\ln \sigma = \ln \sigma_0 - \frac{E_g}{2kT}$$

is with $y = \ln \sigma$ and a linear equation on the type $y = a + bx$, where

$b = -\frac{E_g}{2k}$ is the slope of the straight line.




Fig.2: Conductivity of a semi-conductor as a function of the reciprocal of the temperature

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Hall-effect unit HU 2	11801-01	1
2	Intrinsic conductivity Ge, carrier board	11807-01	1
3	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
4	Digital multimeter, 600V AC/DC, 10A AC/DC, 20 M Ω , 200 μ F, 20 kHz, -20°C... 760°C	07122-00	2
5	Tripod base PHYWE	02002-55	1
6	Right angle clamp expert	02054-00	1
7	Support rod, stainless steel, l = 250 mm, d = 10 mm	02031-00	1
8	Connecting cord, 32 A, 250 mm, red	07360-01	1
9	Connecting cord, 32 A, 250 mm, blue	07360-04	1
10	Connecting cord, 32 A, 500 mm, black	07361-05	2

PHYWE



Setup and Procedure

Setup

PHYWE

The experimental set-up is shown in Fig.1. The test specimen has to be put into the hall-effect-module via the guide-groove. The module is directly connected with the 12 V ~ output of the power unit over the ac-input on the backside of the module.

The plate has to be brought up to the magnet very carefully, so as not to damage the crystal in particular, avoid bending the plate. It has to be in the centre between the pole pieces.

The Hall voltage and the voltage across the sample are measured with a multimeter. Therefore, the sockets on the front-side of the module are used. The current and temperature can be easily read on the integrated display of the module.

The magnetic field has to be measured with the teslameter via a Hall probe, which can be directly put into the groove in the module as shown in Fig. 1. So, you can be sure that the magnetic flux is measured directly on the Ge-sample.

Procedure

PHYWE

At the beginning, set the current I_P to a value of 5 mA. The current I_P remains nearly constant during the measurement, but the voltage changes U_P according to a change in temperature T . Set the display in the temperature mode and be sure, that the display works in the temperature mode during the measurement. Start the measurement by activating the heating coil with the "on/off"-knob on the backside of the module. The specimen will be heated to a maximum temperature of around 145 – 150 °C and the module will stop the heating automatically. Determine the cooling curve of the change in voltage U_P depending on the change in temperature T for a temperature range from 145 °C to room temperature. You will receive a typical curves as shown in Fig. 3 below.





Fig.3: Typical measurement of the probe-voltage as a function of the temperature

PHYWE

Evaluation



Results

PHYWE

With the measured values from Fig. 3, the regression with the expression

$$\ln \sigma = \ln \sigma_0 - \frac{E_g}{2k} \cdot \frac{1}{T}$$

provides the slope $b = (4.05 \pm 0.06) \cdot 10^3 \text{ K}$ (Fig. 4).

With the Boltzmann's constant $k = 8.625 \cdot 10^{-5} \text{ eV}$, we finally obtain

$$E_g = b \cdot 2k = (0.69 \pm 0.01) \text{ eV (Literature value 0.67 eV)}$$




Fig.4: Regression of the conductivity versus the reciprocal of the absolute temperature